## INTENDED USE

### FOR IN VITRO DIAGNOSTIC USE

CLINIQA Liquid QC DipStrip Urinalysis Control is intended for use as an assayed quality control material for various urinalysis reagent strips, qualitative hCG methods as well as confirmatory tests such as Ictotest®.

### SUMMARY

The use of quality control materials to objectively monitor the accuracy and precision of procedures is well established in clinical laboratories. CLINIQA Liquid QC DipStrip Urinalysis Control is provided at two levels to assist in the monitoring of analytical systems within the clinical range.

### PRODUCT DESCRIPTION

CLINIQA Liquid QC DipStrip Urinalysis Control is human urine-based containing constituents of human and animal origin as well as purified chemicals. Preservatives and stabilizers have been added to maintain product integrity. CLINIQA Liquid QC DipStrip Urinalysis Control is a ready-to-use liquid control requiring no reconstitution.

### STORAGE AND STABILITY

CLINIQA Liquid QC DipStrip Urinalysis Control is stable until the expiration date on the tube label when stored unopened at 2-8°C. Do not freeze.

Once opened, tubes of control are stable for 90 days when stored tightly capped at 2-8°C or 20 dipstick immersions have taken place, whichever occurs first. CLINIQA Liquid QC DipStrip Urinalysis Control may be stored for 30 days at 20-25°C.

### PROCEDURE

CLINIQA Liquid QC DipStrip Urinalysis Control should be treated the same as patient specimens and run in accordance with the instructions accompanying the test system (instrument, kit, or reagent) being used. Allow the product to reach room temperature prior to use.

QC materials should be used in accordance with local, state, and/or federal regulations or accreditation requirements.

### INSTRUCTIONS

- Gently invert the tube several times before sampling to ensure homogeneity. This is important to obtain reproducible results.
- For dipstick testing:** remove the cap and dip the reagent strip in the tube as if it were a patient sample. In accordance with the manufacturer's instructions, read the reagent strips either visually or with an instrument reader.

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**For hCG testing:** In accordance with the hCG test kit manufacturer's instructions, use controls as if they were patient specimens.

For the control of confirmatory pregnancy test kits, be sure to use the transfer pipette included in the kit to deliver the correct sample amount to test device.

- Replace the cap immediately after use and store appropriately.

### LIMITATIONS OF PROCEDURE

CLINIQA Liquid QC DipStrip Urinalysis Control should not be used past the expiration date on the tube label. If there is evidence of microbial contamination or excessive turbidity in the product, discard the tube.

CLINIQA Liquid QC DipStrip Urinalysis Control is a stabilized liquid product. To obtain consistent assay values, CLINIQA Liquid QC DipStrip Urinalysis Control requires storage and handling as detailed in STORAGE and STABILITY.

Accurate and reproducible results are dependent upon properly functioning instruments and reagents. The ranges given are intended only as a guideline. Each laboratory should establish their own acceptable ranges and tolerance limits based on their test system.

### EXPECTED VALUES

The expected values printed in this insert were derived from replicate analyses of representative samples of the product and are specific to this lot of CLINIQA Liquid QC DipStrip Urinalysis Control. Testing data used to establish the expected values were derived from multiple instruments. All values have been assigned with manufacturer's reagents available at the time of assay. Subsequent instrument or reagent modifications may invalidate these expected values.

**Refer to [www.cliniqa.com/products](http://www.cliniqa.com/products) for lot specific value assignment. To request a faxed or emailed copy of the value assignment, contact your local distributor or [sma\\_csgroup@bio-technie.com](mailto:sma_csgroup@bio-technie.com).**

### CAUTION

**Human source material. Treat as potentially infectious.**

Each donor unit used in the manufacture of this product has been tested by FDA accepted methods and found non-reactive for the presence of HBsAg and antibody to HIV-1/2, HCV and HIV-1 Ag. While these methods are highly accurate, they do not guarantee that all infected units will be detected. Because no known test method can offer complete assurance the hepatitis B virus, hepatitis C virus, human immunodeficiency virus (HIV) or other infectious agents are absent, all products containing human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

This product contains 0.09% sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form potentially explosive compounds. Flush with excess water upon disposal.

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#### RE-ORDER INFORMATION Liquid QC DipStrip Urinalysis Control Catalog No.

**REF 97201**

Bi-Level, 3 x 2 x 15 mL  
Levels 1 & 2

Made in the USA

Analytes	Siemens Multistix® 10SG Urinalysis Reagent Strips			
	Visual		Siemens Clinitek® Status®/Advantus® Urine Analyzers	
	LOT	Level 1 1804144E	LOT	Level 2 1804145D
Glucose		Negative		100 - 2000 mg/dL or more
Bilirubin		Negative		Small - Large (1+ - 3+)
Ketones		Negative		5 - 160 mg/dL (Tr - Large)
Specific Gravity		1.020 - 1.030		1.000 - 1.020
Blood		Negative		Trace - Large (Tr - 3+)
pH		5.0 - 6.0		6.5 - 8.5
Protein		Negative		Trace - 2000 mg/dL or more (Tr - 4+)
Urobilinogen		Normal (0.2 - 1 mg/dL/EU/dL)		2 - 8 mg/dL
Nitrite		Negative		Positive
Leukocytes		Negative		Trace - Large (Tr - 3+)

Analytes	Siemens Clinitek® Microalbumin Reagent Strips			
	Visual		Siemens Clinitek® Status® Urine Analyzers	
	LOT	Level 1 1804144E	LOT	Level 2 1804145D
Microalbumin		NA		NA
Creatinine		NA		NA

Analytes	Siemens Multistix PRO® 10LS Urinalysis Reagent Strips			
	Visual		Siemens Clinitek® Status®/Advantus® Urine Analyzers	
	LOT	Level 1 1804144E	LOT	Level 2 1804145D
Glucose		Negative		100 - 2000 mg/dL or more
Protein - High		Negative		30 - 2000 mg/dL (1+ - 4+)
Protein - Low		Negative		15 mg/dL (Low)
Ketones		Negative		5 - 160 mg/dL (Tr - Large)
Specific Gravity		1.020 - 1.030		1.000 - 1.020
Creatinine		10 - 200 mg/dL		100 - 300 mg/dL
pH		5.0 - 6.5		6.5 - 8.5
Blood		Negative		Trace - Large (Tr - 3+)
Nitrite		Negative		Positive
Leukocytes		Negative		Trace - Large (Tr - 3+)
Protein		NA		NA
Protein:Creatinine Ratio		NA		NA

Analytes	Roche ChemStrip® Urinalysis Reagent Strips			
	Visual: 10SG, 10UA, 10MD		Roche Urisys 1100®/2400® Urine Analyzers	
	LOT	Level 1 1804144E	LOT	Level 2 1804145D
Specific Gravity		1.015 - 1.030		1.000 - 1.020
pH		5.0 - 6.0		7.0 - 9.0
Leukocytes		Negative		Trace - 2+
Nitrite		Negative		Positive
Protein		Negative		Trace - 500 mg/dL (Tr - 3+)
Glucose		Normal		50 - 1000 mg/dL
Ketones		Negative		Small - Large (1+ - 3+)
Urobilinogen		Normal		1 - 12 mg/dL
Bilirubin		Negative		1+ - 3+
Blood		Negative		Trace - about 250 Ery/µL

Analytes	Roche Combur® Test® Urinalysis Reagent Strips*			
	Visual		Roche Urisys 1100® Urine Analyzer with Roche Combur® Test® UX Strips	
	LOT	Level 1 1804144E	LOT	Level 2 1804145D
Specific Gravity		1.015 - 1.030		1.000 - 1.020
pH		5.0 - 6.0		7.0 - 9.0
Leukocytes		Negative		~10 - 500 Leu/µL (1+ - 3+)
Nitrite		Negative		Positive
Protein		Negative		30 - 500 mg/dL (1+ - 3+)
Glucose		Normal		50 - 1000 mg/dL (1+ - 4+)
Ketones		Negative		10 - 150 mg/dL (1+ - 3+)
Urobilinogen		Normal		1 - 12 mg/dL (1+ - 4+)
Bilirubin		Negative		1+ - 3+
Blood		Negative		5 - 250 Ery/µL (1+ - 4+)

\*Strips not available in the U.S.

Confirmatory and Other Tests			
Test	LOT	Level 1 1804144E	Level 2 1804145D
Ictotest® Reagent Tablets		Negative	Positive
Specific Gravity (refractometer)		1.020 - 1.030	1.000 - 1.020
hCG (Consult*, Genzyme OSOM®, Clinitest®)		Negative (~0 mIU/mL)	Positive (~1000 mIU/mL)
pH paper		4.0 - 6.0	7.0 - 9.0
Total Protein (Sulfosalicylic Acid)		Negative	Positive

\* hCG Dipstick, hCG Urine Cassette, hCG Combo Cassette

32433C\_19 8/24/18

For in vitro diagnostic use

See package insert for proper use

**CLINIQA CORPORATION**


495 Enterprise St.  
San Marcos, CA 92078  
USA

**SALES AND  
TECHNICAL SUPPORT**

P: 800 728 5205  
+1 760 744 1900  
F: +1 760 571 5198  
[www.cliniq.com](http://www.cliniq.com)

**FOR ORDERS AND  
CUSTOMER SERVICE**

P: 800 728 5205  
+1 760 744 1900  
F: +1 760 571 5197  
[sma\\_csgroup@bio-techne.com](mailto:sma_csgroup@bio-techne.com)

CEpartner4U  
Esdoornlaan 13  
3951 DB Maarn, The Netherlands  
[www.cepartner4u.eu](http://www.cepartner4u.eu)

Made in the USA